

IS4000 8mm Base Instruments CDR

September 29, 2011



819005-03 Rev A

UNITED STATES DISTRICT COURT
NORTHERN DISTRICT OF CALIFORNIA

TRIAL EXHIBIT 573-R

Case No. 3:21-cv-03496-AMO

Date Entered

By

Deputy Clerk

Agenda

- Core Team & Design Review Board - *Scott Harrington*
- Market & Product Requirements Overview - *Aaron Carrano*
- Instrument Design
 - Chassis Design - *Ed Donlon*
 - Inputs - *Don Alden*
 - Friction - *Ed Donlon*
- Break! - 11:45am (Get Lunch)
- Clinical Performance - *Maggie Nixon*
- Instrument Performance - *Hsien Hsin Liao*
 - Instrument Hooping performance
 - Comparison with IS3000
 - Modeling of contributors
 - Instrument grip force

Slide 2

ISI Confidential

INTUITIVE
SURGICAL®

Agenda - Continue

- Instrument Performance - continued
 - Safety and reliability - *Niels Smaby*
 - Simulated Life test and sine to death results
 - Wrist back-drive test results
 - Instrument behaviors - *Niels Smaby*
 - Engage trajectories, Grip cal, Cone limits, tool-data from RFID
- *Break! - 15 Minutes*
- V&V Outline - *Craig Tsuji*
- Manufacturing Overview - *Louis Spicciati*
- Current RMA trends - *Marc Bush (Jeff Rampe)*
- Project Risks - *Scott Harrington*
- Schedule Overview - *Scott Harrington*

Slide 3

ISI Confidential

INTUITIVE
SURGICAL®

Instruments Core Team

- **Design Engineering:**
 - Mechanical: Greg Dachs, Ed Donlon, Don Alden, Craig Tsuji, Katie Stoy, Andrew Crews
 - Electrical: Gerry Labonville / Brandon Garbus
 - Systems Analysis: Pushkar Hingwe
 - Software: Brent Tokarchuk
- Project Manager: Scott Harrington
- Product Manager: Aaron Carrano
- Clinical Development: Maggie Nixon
- Customer Service: Jamie Tobin
- Instrument Test: Dallan Schowe
- Manufacturing: Toni Conway
- Manufacturing Engineering: Louis Spicciati / Greg Richmond
- New Product Introduction / Mfg. Test: Garry Doran
- Supplier Development: Hung Pham
- [REDACTED]
- System Test: Steve Auer
- Technical Publications: John Gutierrez
- Purchasing: Chris D'Arienzo / Javed Khan
- Quality: Marc Bush

Slide 4

ISI Confidential

INTUITIVE
SURGICAL®

Instruments Design Review Board

- Clinical Affairs: ~~Victor Chen~~ N/A. No Clinical Trials. (Post-CDR Update)
- Clinical Development: Amy Kerdok
- Customer Service: Brett Robinson
- Mechanical Engineering: Bill Burbank
- Electrical Engineering: Rod Vance
- Manufacturing Engineering (ME): Mike Prindiville
- New Product Introduction: Bob Sundstrom
- NPI Purchasing: David Powers
- Materials: John Hill
- Product Management: Mike Hanuschik
- Product Quality: Krish Krishnanand
- [REDACTED]
- Software Engineering: Greg Toth
- Supplier Development: John Dovala
- Systems Analysis: David Robinson
- System Test: Joe Giallo/Lisa Heaton

Slide 5

ISI Confidential

INTUITIVE
SURGICAL®

IS4000 8mm Base Instruments CDR

Market and Product Requirements



819005-03 Rev A


What's Covered?

Covered	Not Covered
8mm architecture	5mm architecture Advanced Instruments Harmonic ACE Tool Tracking Force Sensing PK Instrument *
	* Working on the PK for Alpha 3.5. This project will have its own CDR.

Slide 7

ISI Confidential

INTUITIVE
SURGICAL®



Instrument Availability

Launch

Strategy

Support Target Specialties



- Hot Shears™ (Monopolar Curved Scissors)
- Permanent Cautery Hook
- Permanent Cautery Spatula
- Maryland Bipolar Forceps
- PK® Dissecting Forceps
- Fenestrated Bipolar Forceps
- Harmonic ACE*
- Curved Bipolar Dissector
- Large Needle Driver
- Mega SutureCut™ Needle Driver
- ProGrasp™ Forceps
- Tenaculum Forceps
- Tip-up Fenestrated Grasper
- Resano Forceps
- Small Graptor™ (Grasping Retractor)
- Large Hem-o-lok® Clip Applier
- Medium Hem-o-lok® Clip Applier
- Atrial Retractor
- Atrial Retractor Short Right
- 5 mm Thoracic Grasper*
- Stapler*
- Vessel Sealer*

**Not covered in this CDR*

Slide 9

ISI Confidential

INTUITIVE
SURGICAL®

Exceptions: full 5mm set.

Instrument Availability

Phase II

Strategy

Complete sets for core



- Mega™ Needle Driver
- Large SutureCut™ Needle Driver
- Black Diamond Micro Forceps
- Cadierre Forceps
- Cobra Grasper
- Double Fenestrated Grasper
- Long Tip Forceps
- Graptor™ (Grasping Retractor)
- Round Tip Scissors

Slide 10

ISI Confidential

INTUITIVE
SURGICAL®

Instrument Availability

Phase III

Strategy

Complete set for cardiac

Mitral
Valve

Revasc

- Micro Bipolar Forceps
- DeBakey Forceps
- Fine Tissue Forceps
- Pericardial Dissector
- Cardiac Probe Grasper
- Potts Scissors
- Small Clip Applier
- EndoWrist Stabilizer
- Snap-fit™ Instrument

Slide 11

ISI Confidential

INTUITIVE
SURGICAL®

Instrument Availability

Retire

Strategy

Discontinue low volume instruments

- PreCise™ Bipolar Forceps
- Thoracic Grasper
- Curved Scissors
- Valve Hook
- Dual Blade Retractor
- EndoPass™ Delivery Instrument

Slide 12

ISI Confidential

INTUITIVE
SURGICAL®

Instrument Availability - Summary*

Launch	22 instruments	Q2 2013
Phase II	10 instruments	Q3 2013
Phase III	9 instruments	Q2 2014
Retire	6 instruments	

**Including advanced instruments, excluding 5mm instruments*


Slide 13

ISI Confidential

INTUITIVE
SURGICAL®

Instrument Types

- Energy Instruments
 - Bipolar Energy
 - Monopolar Energy
 - PK Energy
- Non-energy Instruments



Slide 14

ISI Confidential

INTUITIVE
SURGICAL®

Flush Ports

- Two ports necessary for cleaning
- Co-located with energy connector for Alpha 3
- Strong desire to move flush ports to rear of housing
 - Minimizes overall length when in the ultrasonic cleaner of instrument and simplifies user experience
- Engineering challenge
 - Additional bends in flush tubing reduces pressure



Slide 15

ISI Confidential

INTUITIVE
SURGICAL®

Opening Jaws

- IS3000 instruments have the ability to open jaws of the instrument by manually rotating protruding drive inputs
- IS4000 instrument drive inputs were modified to improve manual rotation



Slide 16

ISI Confidential

INTUITIVE
SURGICAL®

Emergency Release

- Allen key is used to back drive jaws
 - Same as IS3000



Slide 17

ISI Confidential

INTUITIVE
SURGICAL®

Labeling

- Instrument name visible from multiple orientations
 - On robot
 - Off robot
 - Central processing
 - Sterile back table
 - Sterilization tray



Slide 18

ISI Confidential

INTUITIVE
SURGICAL®

End of Life Indicator*

- Utilize extra drive input
 - Rotate drive input to indicate instrument is expired
- Benefits
 - Improve OR efficiency
 - Prevent unnecessary reprocessing



*Users attempt to use an expired instrument **1 in 8** procedures
(over 32,000 occurrences since 01/01/11)*

*Identified as a "should" requirement in PRD

Slide 19

ISI Confidential

INTUITIVE
SURGICAL®


IS4000 8mm Base Instruments CDR

Instrument Design



819005-03 Rev A

Orion Instruments Road Map




Slide 21

ISI Confidential

INTUITIVE
SURGICAL®

8mm 6-cable Architecture Overview

- Substantially the same architecture as presented at PDR
- Changes made to improve performance and manufacturability
- Heavy use of injection molded components
- Utilizes existing IS3000 Distal parts/assemblies




Slide 22

ISI Confidential

INTUITIVE
SURGICAL®

8mm 6-cable Architecture Backend




Slide 23

ISI Confidential

INTUITIVE
SURGICAL®

RFID & Dallas

- Alpha3 Instruments will support both RFID and Dallas,
- Plan of record is to have RFID
- RFID Tag
 - 13mm Square
 - Likely held in with epoxy for A3
 - Will investigate snap or other retaining feature for A4/Beta
- Dallas Board
 - 4 spring loaded pins located on .180" square pattern
 - Held in place with (2) standard Tinnerman washers




Slide 24

ISI Confidential

INTUITIVE
SURGICAL®

Energy Instrument “Disambiguation”




Slide 25

ISI Confidential

INTUITIVE
SURGICAL®

Integrated Energy (Monopolar) Instrument




Slide 26

ISI Confidential

INTUITIVE
SURGICAL®

Integrated Energy (Bipolar) Instrument




Slide 27

ISI Confidential

INTUITIVE
SURGICAL®

Cleaning and Sterilization

- Instrument is designed with standard autoclavable materials
- Instruments will be flushed like IS3000 instruments
- Cleaning verification will start with Alpha 3 Instruments




Slide 28

ISI Confidential

INTUITIVE
SURGICAL®

Magnet (Instrument Presence)

- Use to detect instrument removal
- Instrument will have a magnet similar to IS3000 instruments
- Material is the same, but slightly different geometry (thinner)
- Located under Dallas ID board




Slide 29

ISI Confidential

INTUITIVE
SURGICAL®

Instrument Length




Slide 30

ISI Confidential

INTUITIVE
SURGICAL®

Friction

- Only variable for friction is backend idler pulleys
 - Pulleys need to be as small as possible for packaging, but this adds friction
- Solution: Minimize shaft size and add low-friction coating to pulley and shaft.
- Pulley and shaft geometry not significantly changed from Alpha 2.
- Theoretical worst-case loading could yield 416 shaft, but no problems have been observed.
 - Higher-strength steel available if necessary.
- DLC coated shafts and TiNi pulleys
 - TiNi has same friction characteristics, increased life.
 - Easier to see cables on TiNi pulleys.



Slide 31


ISI Confidential

INTUITIVE
SURGICAL®

Update with photo of TiNi.

Seal Compliance

- Alpha 2 design showed relative motion of shaft and clevis during high grip.
 - “Sewing machine” effect in clinical eval
 - Pitch cables entirely slack under grip
- Alpha 2.5 opened up slot and reduced length of turned-down section to ensure no interference between shaft and clevis seal.
 - Uson test performance remained acceptable
 - Potential for rotation of clevis relative to shaft due to wide slots
- Alpha 3 design decreases length of turned-down section, increases depth of slot.
 - Slot width does not change relative to IS3000



Slide 32

ISI Confidential

Update with photo of TiNi.

Cable Tensions - LND

- Pitch: 5lbs on IS3000 -> 7lbs on IS4000
 - Early tests showed pitch cable slackening during moves under heavy grip
 - Increasing tension prevented this and showed an improvement in performance
 - Pitch cable path does not have reversing bend, so life impacts should be small
- Yaw: 5lbs on IS3000 -> 3lbs on IS4000
 - Keeps overall tension on wrist equivalent during un-gripped moves
 - Does not affect yaw stiffness under heavy grip as grip-open cables are always slack

Slide 33

ISI Confidential


INTUITIVE
SURGICAL®

-Note that pitch performance has a large effect on hooping performance due to the long lever arm between the pitch axis and the tip of the instrument.

-Stiffness gets cut in half when one side of the cable loop goes slack. For pitch this is important as we want the stiffness to be high for optimal tracking. For grip it's less important because when grip is closed hard the grip-open cables are always going to be slack, unless the pretension is insane.

Instrument Inputs

- Orion carriage has continuous rotation.
 - Can use multiple rotations to drive instrument DOFs (Roll, stapler...)
 - *Do not need to control offsets at manufacture*
- Lower input pulley is integral to input disk.
- Upper input pulley is tensioned relative to the lower pulley



Slide 34

ISI Confidential

INTUITIVE
SURGICAL®

Update! (Don Alden) No changes required. DA 26 Sep. Add section view

Input Pulleys

- Keyhole cable termination

- Hypo assemblies arrive with proximal end crimp installed
- Eliminates need for hand-crimper
- Uses same hypo assemblies for all inputs
- Eliminates scrap from trimming tungsten cable after hand crimping.



Slide 35


ISI Confidential

INTUITIVE
SURGICAL®

Update! (Don Alden) No changes required. DA 26 Sep. No changes from PDR. Does not waste Tungsten.

Clamping Pulleys and Shafts

- Bending Stresses:
 - Assume metal post takes all load
 - Larger (0.15" dia) shaft.
 - Assume 49lb load (break strength of cable)
 - Assume plastic adds no strength
 - Maximum von Mises Stress:
 - 68.5ksi → 17-4 PH Stainless Steel (replaces Aluminum)
 - 17-4 PH (H-1150) yield → 126ksi



Slide 36


ISI Confidential

INTUITIVE
SURGICAL®

Update! (Don Alden)

Clamping Pulleys and Shafts

- **Stiffness:**
 - Resolve bending and twisting beam into cable pull
 - Assume stainless steel shaft (0.15" OD)
- 49 Lb causes cable deflection of .0035" if routed to top pulley
- Significantly less (.0016") if routed to lower



Grip close cable always goes to lower input

Slide 37


ISI Confidential

INTUITIVE
SURGICAL®

Update! (Don Alden)

Jaw Cleaning

- Ribs added to jaw inputs to allow opening of jaws for safe cleaning.
- Chassis corners relieved to allow access to jaw inputs.




Slide 38

ISI Confidential

INTUITIVE
SURGICAL®

Instrument Cabling

- Routing adds 1.1" (2.75" total) to tungsten length
 - Compensated by using stiffer cable and using solid rod in place of hypo tube.
 - Reliable manufacturing process for rods has not yet been found.
- Fleet angles for jaw closing improve slightly from 3.7° to approximately 3.0° .



Backend idlers

Slide 39


ISI Confidential

INTUITIVE
SURGICAL®

Update! (Ed Donlon and Katie Stoy) Rods are risky.

Instrument Roll

- Increased RoM
 - Elimination of twist axis allows reduced wrap of rods.
 - Roll ROM of 640° represents an increase over 520° for Gemini.
 - Reducing roll ROM from 640° will improve device life and reduce friction at extreme roll.
- Two-stage hard stop
- Direct drive, 1:1 gearing
 - Simplifies manufacture
 - Decreases part count
- Roll gear alignment marks allow fixed roll offset.



Slide 40

ISI Confidential

INTUITIVE
SURGICAL®

No change but update image Show 640 deg. Tornado? Eliminates burden on roll axis

IS4000 8mm Base Instruments CDR

Clinical Performance



819005-03 Rev A

Alpha2.5 Clinical Performance Testing

Performance tests -

- Instrument Hooping - how much the does the tip deviate when rolling around a fixed point at the tip
- Instrument Power grip - how does the instrument track to bumper position and grasp needles/tissue
- Master Tracking - how well does the instrument track to master inputs (sloppiness)
- Precision suturing - small needles, precise placement and gentle suture interaction
- 'Big bite' suturing - Needle purchase and following with tissue loads, firm suture pulling

Slide 42

ISI Confidential

INTUITIVE
SURGICAL®

Alpha2.5 Testing results

Performance Variable		Alpha2 results (March)		Alpha 2.5 results (September)
Instrument Hooping (LND)	2.5	Straight: 4 – 7mm diameter Bent: n/a	4 3.5	Straight: 1 – 2mm Bent: 4 – 5mm
Instrument power Grip	3	Grip Cal was off Trouble maintaining needle position in jaw		When grip cal working, power grip was equivalent to IS3000
Master tracking (inner axes only)	-		4 3.5	Yaw and pitch Grip (due to delay in grip tracking)
Fine Suturing	3	Some difficulty with tight control of needle tip.	4	Comparable time and stitch precision.
'Big bite' suturing	-	-	4	Needle grip and control was as expected Easily able to follow curve of the needle even with tissue load

1	2	3	4	5
unacceptable	worse	nearly acceptable	equal to IS3000	better than IS3000

Slide 43

ISI Confidential


 INTUITIVE
SURGICAL®

Key Clinical Risks

Category	Detailed Hazard Identification & Rationale	System Behavior & Mitigation	Applicable Functional Area(s)
Imprecise Manipulator motion	Slack and/or friction in instrument (i.e. EndoWrist, snake wrist, single site, etc.) control could result in imprecise motion of instrument wrist causing tissue injury and/or inability to execute necessary surgical tasks.	Instrument movements are uniform (smooth and consistent) over their range of motion. Any backlash, hysteresis, friction, hooping, etc. shall be negligible. Life testing shall be performed to ensure intuitive motion is preserved over life of instrument.	I,M,A
Uncontrolled tip motion	Unintended system executed retraction of the insertion axis of the USM (i.e. due to a supposed instrument removal) could cause tissue injury if the instrument is still present.	System shall redundantly sense presence of instrument to a USM. Retraction of insertion axis for instrument removal only occurs after deliberate, user initiated event (i.e. both instrument release tabs are pressed).	I,M,A,S,E
Non-intuitive	An instrument may break, so that some instrument axes are unable to move as directed by the surgeon. This can cause overall instrument motion to be in an unexpected direction resulting in surgeon confusion and tissue injury.	Instruments shall be designed to withstand typical clinical use and reasonably expected abuse situations throughout the intended life of the instrument. Surgeons should be trained in recognizing the symptoms of unintuitive motion and what can cause it.	I,U

Slide 44

ISI Confidential


 INTUITIVE
SURGICAL®

Key Clinical Risks (con't.)

Detailed Hazard Identification & Rationale	System Behavior & Mitigation	Applicable Functional Area(s)	Detailed Hazard Identification & Rationale
Loss of functionality	Instrument tip (including grips) does not stop or does not stop in a timely manner under a fault condition resulting in tissue injury.	<p>The USM shall stop as soon as possible in accordance to the detection of a system fault. Regardless of USM position or action prior to the fault, the instrument tip shall coast to stop with little or no trauma to any tissue it contacts thus comparable to or less harmful than a dropped laparoscopic instrument.</p> <p>If gripping tissue, a system fault does not release the instrument grips.</p>	M,A,S,E,U
Inability to access patient in emergency conversion	An Instrument gripping /clamping/suctioning or otherwise attached to tissue during a system fault could cause patient injury when it is removed.	The system shall provide a mechanism for quickly (~20 sec) releasing the instrument tips from tissue (i.e. grips, suction, clamping, etc.) during system faults or power losses, so that (a) the instrument releases the tissue, and (b) the sterile adapter releases the instrument.	I,U
Failure of patient contact components	Broken instruments and accessories may generate parts that are difficult to remove, and could cause injury if sharp or become lost in patient anatomy. Potential to result in conversion to open procedure for removal.	<p>Instrument design shall be robust enough for typical clinical use and prevent particles falling into the patient.</p> <p>The User Manual contains appropriate warnings and instructions for safe handling, inspection, set-up & operation.</p>	I,U

Slide 45

ISI Confidential


 INTUITIVE
SURGICAL®

IS4000 8mm Base Instruments CDR

Instrument Performance




819005-03 Rev A

Agenda

- Performance
 - Instrument Hooping performance
 - Comparison with IS3000
 - Modeling of contributors
 - Instrument grip force
 - Instrument tip vibrations
 - Impact of torque ripple
- Safety and reliability
 - Post Fault behavior of instrument tip (HFHF)
 - Simulated Life test and sine to death results
 - Wrist back-drive test results
- Instrument behaviors
 - Engage trajectories, Grip cal, Cone limits, tool-data from RFID

Hooping (video from clinical lab)



- 45 deg. Pitch, 45 deg. Yaw



- 45 deg. Pitch

INTUITIVE
SURGICAL®

Hooping

- Definition...
- Tip Motion is measured by a magnetic tracker.
 - Sensor motion is post-processed to get tool tip motion.




49

ISI Confidential


INTUITIVE
SURGICAL®

mention clinical hooping

Hooping Metric





- Hooping Metric:
 - First project the tip motion onto a plane normal to the initial orientation (wrist pitch at 45 deg.), i.e., plane y-z.
 - Next find the convex hull housing the projected tip motion.
 - The maximum interior distance of convex hull is the metric.



Orion LNDs vs. IS3000 LNDs

- Orion LNDs have similar hooping performance compared to IS3K LNDs.
- Clinical engineers concur during instrument evaluation.
 - However, they usually report 2~2.5 mm LESS error than sensor measurement



51

ISI Confidential


INTUITIVE
SURGICAL®


clinical engineers confirm this with lab sessions.

Source of Hooping Errors

- Want to understand the relative contributions of the following:
 - Measurement Error
 - Uncertainty in magnetic tracker mount position and noise
 - Servo Error
 - Difference between commanded and measured joint positions
 - Carriage and Sterile Adaptor Backlash
 - In a range of 2.5–4.5 degrees (Carriage 1.5–2 deg., Sterile adaptor 1–2.5 deg.)
 - Alpha3 likely has less backlash because ball bearings in the planetary gearbox
 - Actual disc position lags encoder by half backlash
 - Joint Offset
 - Tool “zero” positions not actually zero
 - Wrist pitch can have up to 1.5 deg. offset
 - Instrument Cable Compliance
 - Coupling matrix assumes perfectly stiff instrument cables, which certainly are not
 - Dominant source of hooping error

- A typical symptom: joint hysteresis.
 - Encoder (blue) measures accurate tracking.
 - Magnetic tracker (red) says otherwise.
- A simple cable stretch model can explain this hysteresis...






$$\Delta\theta = \frac{\Delta x}{R} = \frac{T \cdot C}{R} = M \cdot C$$

- $\Delta\theta$ is the disc motion that does NOT contribute to any tool motion.
- Think of an effective disc position that lags the actual position by $\Delta\theta$

- Model estimate the “true” joint position accounting for cable compliance.
- Works very well to explain the hysteresis in both wrist pitch and yaw.




55


ISI Confidential

INTUITIVE
SURGICAL®

Hooping Error



- Run joint angle trajectories through forward kinematics to “simulate” hooping performance.
- With commanded joint trajectories.
 - Max. Interior Distance: 0.032 mm





56

ISI Confidential

INTUITIVE
SURGICAL®

Hooping Error

- With measured joint trajectories.
 - Max. Interior Distance: 0.454 mm





57

ISI Confidential

INTUITIVE
SURGICAL®

Hooping Error

- With measured joint trajectories + backlash model.
 - Max. Interior Distance: 1.194 mm





58

ISI Confidential

INTUITIVE
SURGICAL®

Hooping Error

- With measured joint trajectories + backlash + compliance.
 - Max. Interior Distance: 4.02 mm




59


ISI Confidential

INTUITIVE
SURGICAL®

Hooping Error



- With measured joint trajectories + backlash + compliance + 1.5 deg. pitch offset.
 - Max. Interior Distance: 4.385 mm




60

ISI Confidential

INTUITIVE
SURGICAL®

Hooping Summary

- A typical Orion and Gemini LND at life10:



- Orion hooping performance is comparable to IS3K for LNDs.
- Roughly 1.5 mm un-modeled error remains.

3 bar chart

Grip Force Variation

- Want to maintain a constant grip force.
- Grip force measured by a load cell along with a high-bandwidth amplifier.
- Maneuver:
 - Full grip → yaw to a fixed position → stop.




62

ISI Confidential

INTUITIVE
SURGICAL®

Typical Grip Variation Graphs


- Look at the difference between initial and ultimate grip force.



- Model shows (backup slides) grip variation is equal to the drivetrain friction

Orion vs. IS3000

- 3 new Gemini “blue” LNDs and 6 new Orion A2.0 LNDs.
 - On average, the grip force variation is 40% higher on Orion.
 - Alpha3 carriage friction reduction will reduce the variation.



Grip Variation Summary

- Observe higher grip variation on Orion LNDs.
- Grip variation is exactly the static drivetrain friction during initial grip.
- Consistent with the higher friction measured on Orion.
- Should improve with lower friction Alpha3 design. Dramatic lower friction in the A3.0 motor/gear box.

Instrument Reliability and Safety

- Reliability
 - Simulated life testing
 - Sine cycling to failure
- Safety
 - Back driving the instrument wrist against holding motor torque

INTUITIVE
SURGICAL®

Instrument Reliability, Test Design

- All tests on Production Gemini Instruments and Alpha 2.5 Orion Instruments
- Simulated Life Testing Design(for Large Needle Driver only)
 - Automated Simulated Lives (10 Lives)
 - Includes ultrasonic wash, lube and autoclave for each life.
 - Performance Testing during life span
 - Initial performance (life 0)
 - After first wash/autoclave performance (life 0.5)
 - Then after 1, 5, and 10 lives

INTUITIVE
SURGICAL®

Instrument Reliability, Test Design

- Life Testing Design(for Large Needle Driver only)
 - Simulated Lives
 - 1 Life = From Life Test Protocol for Mega Needle Driver
 - 150 Suture Pulls (more to compensate for no “Tissue Lifts”)



INTUITIVE
SURGICAL®

Instrument Reliability, Test Design

- Life Testing Design(for Large Needle Driver only)
 - Simulated Lives
 - 1 Life = From Life Test Protocol for Mega Needle Driver
 - 150 Suture Pulls (more to compensate for no “Tissue Lifts”)
 - 30 Needle Throws + Pull Throughs



INTUITIVE
SURGICAL®

Instrument Reliability, Test Design

- Life Testing Design(for Large Needle Driver only)
 - Simulated Lives
 - 1 Life = From Life Test Protocol for Mega Needle Driver
 - 150 Suture Pulls (more to compensate for no “Tissue Lifts”)
 - 30 Needle Throws + Pull Throughs
 - 120 Wrist Circles
 - Differences: Instrument changes not performed, No Dips in Pig Blood.



INTUITIVE
SURGICAL®

Instrument Reliability, Test Design

- Life Testing Design(for Large Needle Driver only)
 - Performance Testing during life span
 - Ascension Tracking Measurements
 - Hooping Performance
 - Joint hysteresis
 - Disc Offset Measurement
 - Friction Measurements



INTUITIVE
SURGICAL®


Instrument Reliability, Test Results

- Life Testing Results(for Large Needle Driver only)
 - 6 Orion LND's tested
 - Wrist Pitch Cables Tensioned to 7 lb (vs. 5 lb on Gemini LND)
 - Wrist Yaw Cables Tensioned to 3 lb (vs. 5 lb on Gemini LND)
 - 3 Production LND's tested
 - Failures/Breakages:
 - No Instrument Failures during 10 simulated life test.

INTUITIVE
SURGICAL®

Instrument Reliability, Test Results

- Ascension Tracking Measurements: Hooping Performance
 - Orion on par with Production Gemini




INTUITIVE
SURGICAL®

should stress again that the hooping is good on orion

Instrument Reliability, Test Results


- Ascension Tracking Measurements: Hooping Performance
 - Orion on par with Production Gemini with less variability



INTUITIVE
SURGICAL®

Instrument Reliability, Test Results

- Disk Offset Changes: Group Comparison, Disks 3 and 4 (yaw/grip)
 - Consistently more movement on Orion Disk offsets. (w/ one outlier)



■ Grip Cal Compensation should result in no detrimental performance effects.

- Hooping not affected due to symmetry of drift.

INTUITIVE
SURGICAL®

Instrument Reliability, Test Design

- Grip Max and Sine Cycle Yaw until distal cable breakage
 - To test common failure mode at distal idler pulleys.

INTUITIVE
SURGICAL®

Instrument Reliability, Test Results

- Grip Max and Sine Cycle Yaw (discrete cosine trajectory start/stop, check) until distal cable breakage

Gemini Failures	Orion Failures
641 cycles	660 cycles
605 cycles	457 cycles
579 cycles	436 cycles

Note: The grip torque at the motor for Orion instrument was at 0.18 Nm compared to 0.16 Nm for IS3000. However, the grip force generated was comparable.

INTUITIVE
SURGICAL®

Instrument Reliability, Results

- Why do the Orion Instruments break sooner in Max Grip Sine Cycle?
 - Orion Stiffer Drive Train May cause higher wrist yaw peak torques?
 - Gemini back end friction increase during test “protecting” wrist?
 - **RISK: Priority for further investigation.**
 - Measure grip with grip meter during test.
 - Determine fastest acceleration on WY in following
 - Match Sine or Cos Trajectory to be within reasonable acceleration bounds
 - Tune maximum grip torque on Orion to better match IS3000 grip forces.
 - Sensitivity study (to max grip torque).

Instrument Safety, Back driving Pitch and Yaw

	Pitch back driving result	Yaw back driving result
IS3000	Reached ROM limit but did not break against holding torque of 0.4 Nm at the motor.	Yaw cable broke at the wrist beyond 0.35 Nm motor holding torque (tested in increments of 0.05 Nm holding torque at the motor)
Orion	Reached ROM limit but did not break against holding torque of 0.45 Nm at the motor.	Yaw cable broke at the wrist beyond 0.35 Nm motor holding torque



Result of back driving the pitch



Result of back driving the yaw


INTUITIVE
SURGICAL®

Instrument Behaviors and Algorithms

- Instrument engagement behaviors
- Grip Cal
- Cone Limits
- Instrument Data, Disc Offsets, Plug-N-Play

INTUITIVE
SURGICAL®


Carriage and Instrument Behaviors and Algorithms, Instrument Engagement Behaviors



Arbitrary disk offsets on instruments: unable to ensure that we start "on one side" of engagement.
 wrist pitch, wrist yaw and grip engage only when roll is stopped or moving slowly (cannula scraping). (+0.5 sec)
 roll, an additional full rotation on roll(+ 0.25 sec), followed by second behavior to go to zero roll (+0.5 sec) (Eliminated by fixed roll disk offset)
 120 deg greater ROM on roll. (+0.08 sec)

Carriage and Instrument Behaviors and Algorithms, Grip Cal

- Finds the grip close position
- “Torque Ripple” in Alpha 2 Carriages prevented accurate determination of tool bumper. **Tuning deferred to Alpha 3.**



INTUITIVE
SURGICAL®

Carriage and Instrument Behaviors and Algorithms, Wrist Cone Limits

- To ensure intuitive virtual slave limit feedback to masters, Wrist Cone Limits are now turned on for all Orion Instruments (IS3000 only snake wristed, LND & BD Micro)
 - Results in reduction in wrist range of motion (inside the joint limits, especially when close to the remote center)
 - Example: Monopolar Spatula instrument
 - Joint signal ranges: WP: 80 deg, WY: 85 deg.
 - At 0.25m insertion, Wrist Cone Limit = ~78 degrees
 - At 0.035m insertion, Wrist Cone Limit = ~65 degrees.
 - Loss of 20 deg of WY motion = loss of ~6mm of reach back
 - Algorithmic fix needed.
 - IS3000 does work well.

Carriage and Instrument Behaviors and Algorithms, Instrument Data/RFID Tags

- Header Data (programmed into RFID and tested to work)
 - Tool UID
 - Laser ID/Tool SN (used for energy mapping)
 - Disc Offsets
- Plug-N-Play Data (not yet programmed and tested)
 - Plenty of memory space on the RFID tags.
 - Not believed to be a risk

INTUITIVE
SURGICAL®

Orion Instrument RFID Tags

Testing Summary

- Two different RFID tag companies are developing a RFID tag for ISI
- The next round of autoclave testing will be done in the CT Office and will include cleaning and sterilization cycles that are done in a Pre-Vacuum autoclave and that are not batched
 - This testing will be the most similar to what conditions the parts will see in the field
 - Testing with a read/write every 5 cycles at a nominal read distance
- Each company can provide us with a tag(s) that can fulfill different requirements for Orion

Company	Protocol: ISO 15693	Protocol: ISO 14443-B	Preferred Atmel Chip	Autoclaveable	Gamma Stable	Application: Instrument vs. Cannula
RCD Technology		X	X	X		Instrument
NewAge Industries	X			X	X	Both
NewAge Industries		X	X	X		Instrument

- Future Testing Plan:
 - 1000 samples received from RCD and are ready to start testing around 10/3/11
 - 1000 samples due from NewAge (ISO 15693) on 11/11/11
 - 1000 samples from NewAge (ISO 14443-B) with preferred Atmel chip is still TBD because of issues for NewAge sourcing chip from Atmel

Slide 85
ISI Confidential

INTUITIVE
SURGICAL®

IS4000 8mm Base Instruments CDR

Verification and Validation Outline



819005-03 Rev A

Top Level Test Plan

REQUIRED V/V TESTING FOR IS4000 INSTRUMENTS:

- **Design Verification** - *One Protocol covering general instrument architecture plus individual protocols for each instrument.*
- **Design/Clinical Validation** - *One Protocol covering general instrument architecture plus individual protocols for each instrument.*
- **Software Compatibility/Validation**
- **Reliability/Life Testing** - *Separated into 6 protocols, based on intended function of the instrument. (1) Monopolar Instruments, (2) Bipolar Instruments, (3) Needle Driver Instruments, (4) Grasper Instruments, (5) Clip Applier Instruments, and (6) Specialty Instruments*
- **Electrical Safety and EMC Testing**
- **Cleaning Validation**
- **Sterilization Validation**
- **Biocompatibility** - *No new materials at this time.*
- **Surgeon Evaluation**
- **Packaging Validation** - *Longer instrument will require new box*
- **Human Factors**

Slide 87

ISI Confidential

INTUITIVE
SURGICAL®

Architectural Requirements 823033-40 (same format as IS3000 ARD)

Major changes from IS3000 ARD:

- Update Reference section
- Instrument over-all dimensional changes:
 - Width of backend narrower
 - Height of backend larger
 - Length 2" functionally longer past remote center (except for Harmonic Ace)
 - New foot print for Instrument-to-ISA disks
- Adding fifth input
- Requirements for RFID function added
- Established new "datum" to measure length of instrument
 - Was from centerline of distal input disks to distal tip of instrument
 - Now is from bottom of chassis to distal tip of instrument
- Use 5th axis as a end of life indicator

Slide 88

ISI Confidential

INTUITIVE
SURGICAL®

Update! Craig Tsuji

Functional Requirements 823123-40 (same format as IS3000 FRD)

Major changes from IS3000 FRD:

- Increased Roll ROM
- Removed the Clinical Requirements
- Otherwise, document largely remain the same as IS3000 FRD

Slide 89

ISI Confidential

INTUITIVE
SURGICAL®

Scott H.

IS4000 8mm Base Instruments CDR

Manufacturing Overview



819005-03 Rev A

Part Numbers

- Planning BOM PB-OR-A3-INST
 - Distal Assemblies (same as current instruments)
 - 18 instruments
- Planning BOM PB-INST-A3NEW
 - Proximal Assembly (common parts)
 - 19 parts
- Assy BOM
 - 7 Instruments at -01X

Number	Description
470006-01X	ASSY,LARGE NEEDLE DRIVER,8MM,IS4000,ENG
470033-01X	ASSY,BLACK DIAMOND MICRO FORCEPS,8MM,IS4000,ENG
470172-01X	ASSY,MARYLAND BIPOlar FORCEPS,8MM,IS4000,ENG
470179-01X	ASSY,MONOPOLAR CURVED SCISSORS,8MM,IS4000,ENG
470183-01X	ASSY,PERMANENT CAUTERY HOOK,8MM,IS4000,ENG
470189-01X	ASSY,DOUBLE FENESTRATED,8MM,IS4000,ENG
470205-01X	ASSY,FENESTRATED BIPOlar FORCEPS,8MM,IS4000,ENG

Slide 91

ISI Confidential


 INTUITIVE
SURGICAL®

Manufacturing Strategy

Instrument Front-end VS Instrument Backend

- Leverage current processes where available at front-end:
 - Common grips IS2000 provide volume & cost benefits.
 - Production methods benefit from known, mature processes.
 - Similarities to IS2000 reduce risks at launch.
- Focus on changes at Backend
 - DFM advantages over IS2000
 - No handcrimping of cables, auto tensioning, no re-tensioning required
 - Requires development of new assembly fixtures.
 - RFID requires new “DCP” like programming solution.
 - IPT Tester - new fixture to support 5 axis and process improvements.

Slide 92


ISI Confidential

INTUITIVE
SURGICAL®

Luis S and Greg R

Process Development / Assembly Tooling

IS2000 STANDARD PROCESS




Slide 93

ISI Confidential

INTUITIVE
SURGICAL®

Process Development / Test Process




Slide 94

ISI Confidential

INTUITIVE
SURGICAL®

Tooling Requirements

- **Assembly Tooling**
 - **Cable Combs** - Same at IS2000 but longer
- **MINI-Brick** -
 - Hold / Position Instrument
 - Restrain Input rotation
- **Bearing presses** - install bearings
- **Tensioner** - set initial Tension
- **Snap ring tool**



Slide 95

ISI Confidential

INTUITIVE
SURGICAL®

Tooling Requirements - cont

■ RFID Programmer

- New HW / SW System
- Hardware to interface to instrument / RFID
- SW to support required functionality
- Functionality of current “DCP” -
 - Program and verify instruments
 - Read instrument information for IPT, Marker, and Labeling system.
 - Write and verify process control flags.
 - Maintain security of programming system
- NEW - Write instrument OFFSETs to RFID Chip
- NEW - Read Energy Dallas ID and write to RFID Chip

■ Resources

- ME - Mechanical layout
- TE - HW
- TE - Diagnostics (SW)

Slide 96

ISI Confidential

INTUITIVE
SURGICAL®

Orion IPT

- New HW / SW System
- 5- Axis drive capability
- Functionality of current IS1200/2000 IPT -
 - Read instrument and load test scripts
 - Offset Measurement
 - Range of Motion/Friction
 - Cabling Verification
- NEW - Write OFFSETs to RFID Chip or DCP
- NEW - Automatic Cable Tensioning
 - Consistent cable tensioning
 - Allows use of advanced techniques
- Next Phase is S/W integration by mid/late Oct

- Resources
 - ME - Mechanical layout
 - TE - HW / Controls
 - Test Eng Diagnostics (SW)



Slide 97

ISI Confidential

INTUITIVE
SURGICAL®

Production Line Lifecycle Development

Alpha

- Alpha 3 (Q4) - IPT and DCP offline - develop Orion specific line
- Alpha 4 (Q1) - IPT and DCP online - initiate EFT qualification

Beta

- Complete final EFT qual
- Perform process validation

Gamma

- Duplicate Orion line(s) in support of launch volumes
- Launch builds in Sunnyvale

Production

- Demonstrate stability
- Transfer to Mexicali plant -

Slide 98

ISI Confidential

INTUITIVE
SURGICAL®

Supply Base

- Prime components follow current IS2000/3000
 - Cables / Hypotubes - Baird Industries
 - Grips / Clevis - Britt / Indo-MIM / Remmele
 - Maintubes - Easton
 - Injection Molding - UPG

- New Parts/Suppliers
 - Hypotubes - New process at Baird
 - RFID - New process at RCD and New Age Industries
 - Clamping Pulley - Swiss Precision

Slide 99

ISI Confidential

INTUITIVE
SURGICAL®

Process Validation Plans

- Process Validation to follow standard instrument format
 - Demonstrate EFT functionality and PMC
 - Demonstrate availability of work instructions and trained technicians.
 - Execute three build lots (10 units minimum) for each instrument type.
 - Evaluate In-process Testing results.
 - Gather Variable Data and Assess margin to spec
 - Evaluate product quality against specification. (PQA functional test)
 - Evaluate Labeling and DHR documentation.

Slide 100

ISI Confidential

INTUITIVE
SURGICAL®

Material Costs

MATERIAL COST			
	IS2000	IS4000	Com
LND	\$ 157	\$ 151	\$ 89
Bipolar	\$ 181	\$ 180	\$ 128
MCS*	\$ 270	\$ 269	\$ 135

* MCS includes cost of 1 Tip Cover

- Instrument COGS is driven by Grip and Cable costs.
- Hypotube cost estimate may increase due to process limitations
- IS4000 cost reduction is sensitive to cost of bearings.
- Moderate improvement to standard Labor & O/H by DFM.
- Stored Offsets to have significant impact on Mfg efficiency (first pass yield)

Slide 101

ISI Confidential

INTUITIVE
SURGICAL®

IS4000 8mm Base Instruments CDR

Current RMA Trends



819005-03 Rev A

Predecessor Instrument RMA


Instrument FA
(Normalized by Procedure)

> 0% > .10% > .20% > .30% > .40% > .50%

Complaint YES

Failure Analysis	Q3-2010		Q4-2010		Q1-2011		Q2-2011		Trend	Notes	Action Tracking	Owner
	Proc.	QTY	Proc.	QTY	Proc.	QTY	Proc.	QTY				
TOTAL		1398	2.01%	1301	1.67%	1343	1.65%	1402	1.59%	—		
NTF (No Trouble Found)	341	0.49%	339	0.43%	390	0.48%	367	0.41%	↓	PT under development for iRMA team	N/A	GirishS
Cables (Broken/Frayed)	289	0.42%	296	0.38%	299	0.37%	293	0.33%	↓			
Scissors (Cutting Performance)	172	0.25%	138	0.18%	134	0.16%	162	0.18%	—			
Grips (Broken/Bent)	121	0.17%	102	0.13%	133	0.16%	153	0.17%	—			
Housing (Broken)	124	0.18%	112	0.14%	126	0.15%	110	0.12%	↓	Design update (materials) implemented	C45185	SharkS
Energy Pins/Plugs (Damaged)	41	0.06%	40	0.05%	47	0.06%	64	0.07%	↑	Design update (loc/ite/chassis) in-progress	C52176	LewisI
Bipolar Conductor Wire (Broken)	35	0.05%	41	0.05%	27	0.03%	46	0.05%	—	Supplier update (cables) under investigation	C53212	MikeB
Tube Extension (Broken)	49	0.07%	34	0.04%	45	0.06%	37	0.04%	↓	Design update (3 to 4 keys) implemented	C43819	MikeB
Harmonic Inserts (Ace/Shears)	25	0.04%	32	0.04%	14	0.02%	24	0.03%	—	Shears (Clamp Arm) Ace (Broke Blade/Self-test)		
Main Tube (Scratches)	45	0.06%	46	0.06%	29	0.04%	21	0.02%	↓			
Luer Plate (Broken)	36	0.05%	37	0.05%	28	0.03%	19	0.02%	↓	Design update (-05) implemented	C46402	SharkS
Clamping Pulleys (Loose/Cracked)	9	0.01%	8	0.01%	8	0.01%	17	0.02%	↑	Torque specification added to drawings	C52702	SharkS
Hooks & Spats (Broken)	18	0.03%	22	0.03%	15	0.02%	15	0.02%	↓	Design update (-06) implemented	CAPA1339	MikeB
Other	93	0.13%	54	0.07%	48	0.06%	74	0.08%	—			

Overall complaint rate stable (flat) - several design improvements shipping to customer base.
Harmonic inserts monitored for Shears vs Ace failure modes and rates.



Slide 103

QRB 2011 Q2

INTUITIVE
SURGICAL®

Orion Instrument Comparison

Historical Failures	Orion Mitigations
No Trouble Found	No specific mitigations planned
Drive Cable failure	Improved Cable tension control
Scissor cutting performance	Cut Test Validation
Grips Bent/Broken	Life Test Validation
Housing Broken	Using Radel material
Energy Plugs Damaged	New connector design minimizes the pin exposure
Conductor Wire Break	Life Test Validation
Tube extension broken	Life Test Validation
Harmonic Inserts	N/A for Orion standard instruments
Main Tube scratches	Easton Tube material
Luer Plate broken	Cleaning and Life Test Validations
Clamping Pulleys Loose/cracked	Life Test Validation
Hooks-spats (Plastic clevis)	Architectural Requirements for axial load Non-glass-filled Ultem 1000 distal clevis

Slide 104

ISI Confidential


 INTUITIVE
SURGICAL®

Highlight the effective specs. Drive cable failure and Energy Plugs Damaged

IS4000 8mm Base Instruments CDR

Project Risks



819005-03 Rev A

Project (Technical/Schedule) Risks

Item	Mitigation/Update	Risk
Excess hooping due to compliance and/or backlash	- See Performance Slide	Low
Long instrument - increased potential for capacitive leakage and deflection	Waiting for Alpha 3 build to test this issue Use shorter length	Low
Clinical usage of Harmonic Ace (shorter compared to other IS4000 instruments)	Work with Ethicon to develop longer Harmonic Blade Insert Use shorter length	High (if we absolutely need a longer length Harmonic Blade Insert)
Increase Roll Angle: Life Impact	Need to evaluate	Medium

Slide 106

ISI Confidential


 INTUITIVE
SURGICAL®

Project (Technical/Schedule) Risks (con't.)

Item	Mitigation/Update	Risk
Increase engagement time	Reducing engagement time to 1.2 seconds by eliminating roll disc offset	Low
Earlier Failure from sine cycling	More testing/ investigate the root cause	High
RFID	Sterilization/life testing Reliability testing Dallas chip as a backup	Medium
Hypo-tube manufacturing	Alternate designs are being considered Improve manufacturing process Use IS3000 Hypo-tube	Medium

Slide 107

ISI Confidential


 INTUITIVE
SURGICAL®

IS4000 8mm Base Instruments CDR

Project Schedule



819005-03 Rev A

Orion Instrument Schedule

2011																																																			
18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52																	
M	M	M	M	M	J	J	J	J	J	J	J	J	A	A	A	A	A	S	S	S	S	O	O	O	O	O	N	N	N	N	D	D	D	D																	
2	9	16	23	30	6	13	20	27	4	11	18	25	1	8	15	22	29	5	12	19	26	3	10	17	24	31	7	14	21	28	5	12	19	26																	
Alpha 3 design																					Fabrication											g	assem	3	FVI test																
																					Alpha 4														11																

- Ordered long lead items - still waiting for quotes
 - Ordered machined assembly for end of October delivery (Qty: 12)
- Total ordered quantity: 500 Instruments
 - Quantity needed for testing: 231
 - Initial build quantity: Approx. 100
 - Remaining will be built per demand
- Mini Brick: due on third week of Oct.
- IPT (Phase 1): due on second week of Nov.
- Alpha 4 design will start third week of October
 - Concepts for the end of life indicator (Fifth Input)
 - Industrial Design
 - Harmonic Ace
- Technical reviews: End of November and December
- Alpha 4 design release: End of December

Slide 109


ISI Confidential


 INTUITIVE
SURGICAL®

4/12/2021

Scott Harrington

Orion Instrument Schedule



- Alpha 4 fabrication: first week of January, 2012
- Finalize feasibility assessment of, fifth input (end of life indicator), and moving of the flush ports (back of the instrument).

■ [REDACTED]

■ [REDACTED]

■ [REDACTED]

■ [REDACTED]

- Finalize Beta Design: End of April, 2012
 - Technical reviews: End of March and April
- Finalize V&V strategy: End of February, 2012

Slide 110


ISI Confidential

 INTUITIVE
SURGICAL®

4/12/2021

Scott Harrington

Orion Instrument Schedule



- [REDACTED]
- FDR: Third week of January
- Phase 2 Instrument Development: Starts third week of January
- [REDACTED]

Slide 111

ISI Confidential

 INTUITIVE
SURGICAL®

4/12/2021

Scott Harrington

INTUITIVE
SURGICAL®

Trial Ex. No. 573-R, Pg. 112 of 124

Instrument Action Items from PDR - Resolved

No.	Issue/Action	Resp.	Target Completion	Final resolution	Closed
1	Evaluate optimal direction of cautery cord as it comes off the instrument housing	Greg Dachs	CDR	Cautery connector board incorporates the optimal position. CDE provided recommendation after meeting with team to review Projet prototype.	Oct. 2011
2	Add 10 - 5mm reducer to Accessory launch deliverables	Joe Orban	CDR	PDR combined Instruments and Accessories. There will be a separate CDR for Accessories and this item will be addressed at that CDR.	N/A for Instruments

(Post-CDR Update)

Slide 113

ISI Confidential

INTUITIVE
SURGICAL®

Questions?

114

ISI Confidential

INTUITIVE
SURGICAL®

Thank You For Your Time!
Please sign the attendance
sheet

Slide 115

ISI Confidential

INTUITIVE
SURGICAL®


Backup Slides

116


ISI Confidential

INTUITIVE
SURGICAL®


Multi-Dof Model Validation



- Simultaneously exercising wrist roll, pitch and yaw.
 - Model gives estimated joint pos., which are run through forward kinematics.
 - Project tip motion on to the “y-z plane”.



- A simple mass-spring model can explain the grip variation in “grip and step yaw”




- Reducing friction helps reduce grip force variation seen in “grip and step yaw”.

118


ISI Confidential

INTUITIVE
SURGICAL®

dropped



- The simple mass-spring model also says the following:
 - It is important to have balanced friction between the 2 yaw jaws.
 - Easier to achieve on Orion with disc offsets
 - It is important to maintain a constant grip torque
 - Remaining torque ripple is why Orion shows larger grip force variation during the “grip and sine yaw” move
 - “Spring” length relates to spring force
 - The grip dof position tells us a lot about the grip force variation.




120

ISI Confidential

INTUITIVE
SURGICAL®

dropped


Instrument Reliability, Backdriving the pitch and the yaw



ISI Confidential

INTUITIVE
SURGICAL®

Instrument Reliability, Backdriving the pitch and the yaw



ISI Confidential

INTUITIVE
SURGICAL®

Instrument Bringup Status

- Significant additional Software/Parameter effort still required.

EndoWrist Monopolar Cautery Instruments				WDEF CREATED	PARAMETERS TUNED	SW/TUNING EFFORT
420179	Hot Shears™ Monopolar Curved Scissors		Launch	YES	NO	Large
420183	Permanent Cautery Hook		Launch	YES	NO	Small
420184	Permanent Cautery Spatula		Phase 2	YES	NO	Small
EndoWrist Bipolar Cautery Instruments						
420172	Maryland Bipolar Forceps		Launch	YES	NO	Large
420227	PK® Dissecting Forceps		Launch	NO	NO	Small
420205	Fenestrated Bipolar Forceps		Launch	YES	NO	Small
420344	Curved Bipolar Dissector		Launch	NO	NO	Large
Ultrasonic Energy Instruments						
420275	Harmonic ACE® Curved Shears 8 mm		Launch	NO	NO	Small
EndoWrist Needle Drivers						
420006	Large Needle Driver		Launch	YES	YES	Medium
420194	Mega™ Needle Driver		Launch	NO	NO	Small
420309	Mega SutureCut™ Needle Driver		Launch	NO	NO	Medium


EndoWrist Graspers				MOEF CREATED	PARAMETERS TUNED	SW/TUNING EFFORT
420033	Black Diamond Micro Forceps		Phase 2	YES	NO	Small
420093	ProGrasp™ Forceps		Launch	NO	NO	Large
420049	Ladriere Forceps		Phase 2	YES	NO	Small
420207	Feneculum Forceps		Launch	NO	NO	Small
420036	DeBakey Forceps		Launch	NO	NO	Small
420189	Double Fenestrated Grasper		Phase 2	YES	NO	Small
420600	Tip-up Fenestrated Grasper		Launch	NO	NO	Small
420318	Small Graptor™ (Grasping Retractor)		Launch	NO	NO	Medium
Clip Applier						
420230	Large Hem-o-lok® Clip Applier		Launch	YES	NO	Large
420327	Medium Hem-o-lok® Clip Applier		Launch	NO	NO	Large
Specialty Instruments						
420204	Atrial Retractor		Launch	NO	NO	Medium
420246	Atrial Retractor Short Right		Launch	NO	NO	Medium
420343	5 mm Thoracic Grasper		Launch	NO	NO	Large

INTUITIVE
SURGICAL®

Orion Instrument RFID Tags

Testing Summary

- 6 different RFID tags from 3 different Manufacturers were tested to 301 cleaning cycles
 - 1 cleaning cycle = 15min ultrasonic cleaner, 1 autoclave cycles (no preVacuum) at 132°C for 4min (testing completed at test lab) batched into groups of 5



Tag Mfgr & Type	# of test samples	Failures @ cycle #
<u>RCD Technology</u> † (Handmade HF 5010 Rev1 tag)	11	56, 196, 224, 245, 273, 273, 273*
<u>Syama Tech</u> (7mm Epoxy Autoclave)	10	140, 140, 154, 161, 182, 182, 196, 196, 196, 196
<u>NewAge Industries</u> † (RFID-HT-16mm High Temp Tag)	10	Zero Failures!
<u>NewAge Industries</u> (Part# NS1 10mm)	10	91, 161, 161, 161, 189, 280, 301
<u>NewAge Industries</u> (CA7991026 -7.5mm PCB Tag)	3	77, 105, 147
<u>Atmel</u> PN:174305	10	21, 49, 49, 49, 56, 70, 70, 84, 98

Slide 124
ISI Confidential

*Failure analysis showed failures were unique to chip used in prototypes

Continued to next round of prototyping

INTUITIVE
SURGICAL®