

# **Exhibit 1**



# da Vinci® EndoWrist® Repairs

## **Stop throwing away money on your da Vinci® EndoWrist® devices!**

SIS can now service your da Vinci® devices including repair and use counter reset.

Important facts:

- The da Vinci® EndoWrist® is a “multi-use” medical device. Multi-use devices, such as endoscopic instruments, have always been eligible for repair.
- The repair of da Vinci® EndoWrist® does not alter the intended use, method of use, functionality or performance of the device in any way.
- The da Vinci® Robot interacts with the repaired EndoWrist® identically and the robot cannot be affected by the repaired device in any way.
- The robot communicates with the EndoWrist® prior to surgery only. Prior to surgery, the robot confirms the serial number, model number and remaining uses. The repaired device will function identically to the new OEM EndoWrist®.
- A repaired EndoWrist® is not an alternative or replacement device. It is an original da Vinci® manufactured device that has been repaired to original specifications.

## **It's time to challenge the status quo**

SIS da Vinci® EndoWrist® repair services will boost your bottom line and help provide the prompt and tailored responses your OR demands.

Call your local SIS representative or our Corporate Office at 800.747.8044.